

INFINITY

COLLECTION

TRUE-INFINITE POSITION SEATING

The Ultimate in Relaxation and Support

The Infinity Collection from Pride® Power Lift Recliners features true-infinite position seating and pocketed coils designed to deliver comfort, style and performance. Available in four different fabric selections, the Infinity Collection models deliver maximum relaxation and will match any décor.



USB Charger Port (standard)



Scan here for more information on the **Infinity Collection**®



FOLLOW US ON SOCIAL MEDIA

PRIDE®
MOBILITY
LIVE YOUR BEST®

FEATURES & BENEFITS

- Pocketed coils deliver individualized support and relaxation
- Premium furniture grade laminate/hardwood frames for added strength and durability
- Easy-to-use hand control with USB charger (not available with optional heat and massage)
- Integrated battery backup ensures chair is operational when power goes out
- Energy-efficient power box equipped with self-tested electronics and 9' cord
- Standard lumbar pillow
- Large, dual pockets for extra storage
- Quiet and smooth lift system

OPTIONS*

- Deluxe heat and massage
- Residential cleaning kit
- Head pillow
- Protective pad

* Additional charges apply for all options. Extra lead time may apply for some options.

MODEL	LC-525iPW	LC-525iS	LC-525iM	LC-525iL
POSITION	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'2" & Below	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'3" & Below	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'4" - 6'0"	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'9" - 6'2"
SIZE	Petite Wide	Small	Medium	Large
WEIGHT CAPACITY	375 lbs.	375 lbs.	375 lbs.	375 lbs.
SEAT-TO-FLOOR	19"	19"	20"	21"
SEAT DEPTH	19"	19"	20"	22"
SEAT WIDTH	24"	19"	20"	22"
TOP OF BACK TO SEAT	28"	26"	28"	30"
BACK STYLE	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes
DELUXE HEAT AND MASSAGE	Optional	Optional	Optional	Optional
SEATING TYPE	Pocketed Coil Springs	Pocketed Coil Springs	Pocketed Coil Springs	Pocketed Coil Springs
DISTANCE FROM WALL	19"	20"	20"	21"
OVERALL WIDTH	38"	34"	34.5"	37"
CHAIR WEIGHT	130 lbs.	118 lbs.	130 lbs.	130 lbs.
WARRANTY	7-Year Pro-rated ¹	7-Year Pro-rated ¹	7-Year Pro-rated ¹	7-Year Pro-rated ¹



(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1. Lifetime warranty on steel lift frame.

3. Quick Ship fabric

2. Five-year warranty on moisture barrier (void if punctured).

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]