

|  |   |
|--|---|
| <b>Model</b>   | Go-Go® Ultra X  |
| <b>Part Number</b>                                   | 3-wheel: SC40X<br>4-wheel: SC44X  |
| <b>Weight Capacity</b>                               | 260 lbs.  |
| <b>Maximum Speed <sup>1</sup></b>                    | Up to 4 mph   |
| <b>Ground Clearance <sup>3</sup></b>                 | 3-wheel: 3.1" at deck<br>4-wheel: 3" at deck  |
| <b>Turning Radius <sup>3</sup></b>                   | 3-wheel: 33" 4-wheel: 44"   |
| <b>Overall Length <sup>3</sup></b>                   | 3-wheel: 37" 4-wheel: 39.5"   |
| <b>Overall Width <sup>3</sup></b>                    | 19.5"   |
| <b>Seat-to-Deck Height</b>                           | 15.75" - 17.25"   |
| <b>Seat-to-Ground Height</b>                         | 18.25" - 20.25"   |
| <b>Front Tires</b>                                   | 3-wheel: 2" x 8" solid<br>4-wheel: 2" x 7" solid  |
| <b>Rear Tires</b>                                    | 2.5" x 8" solid   |
| <b>Range Per Charge <sup>1,2</sup></b>               | 3-wheel: Up to 6.9 miles<br>4-wheel: Up to 7.2 miles  |
| <b>Total Weight Without Batteries</b>                | 3-wheel: 76 lbs.<br>4-wheel: 80 lbs.  |
| <b>Heaviest Piece When Disassembled <sup>3</sup></b> | 3-wheel: 29 lbs. (front section)<br>4-wheel: 34 lbs. (front section)                                      |
| <b>Standard Seating</b>                              | Type: Compact, foldable<br>Weight: 23 lbs.<br>Material: Black vinyl<br>Dimensions: Width: 17", Depth: 17" |
| <b>Motor Output</b>                                  | 24-volt DC motor,<br>sealed transaxle   |
| <b>Brakes</b>  | Regenerative and<br>electromechanical   |
| <b>Battery Requirements <sup>4</sup></b>             | (Two) 12 AH<br>Weight: 9 lbs. each  |
| <b>Battery Charger</b>                               | Off-board, 2 amp  |

## HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY  
DME PDAC OR DOES NOT MEET CRITERIA

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup> Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>4</sup> Battery weight may vary based on battery manufacturer.

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

GOGOULTRAX\_061819

**GOGO**<sup>®</sup>  
**ULTRAX**

**ECONOMY  
CLASS PRICE...**

**...FIRST CLASS  
FEATURES!**



**LIVE YOUR BEST<sup>®</sup>**

pridemobility.com |    

**ECONOMY CLASS PRICE...**  
**...FIRST CLASS FEATURES!**



The Go-Go® Ultra X is loaded with features such as one-hand disassembly and a convenient drop-in battery box for worry-free travel. With a 260 lb. weight capacity, a maximum speed up to 4 mph, a per-charge range up to 6.9 miles on the 3-wheel and 7.2 miles on the 4-wheel, Go-Go Ultra X is an exceptional value.

**Additional Features**

- Convenient dual voltage charger can charge battery pack on-board or off-board
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability
- Pride's exclusive black, non-scutting tires

**Color**




Red

**Accessories\***

- Cane/crutch holder (single, double)
- Cup holder
- Double crutch holder
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Battery pack

**Military Patches**



Convenient drop-in battery box for worry-free travel.



Feather-touch disassembly for quick transport or storage.



Personalize your Go-Go® with one or more of our optional accessories.

\* Only one rear accessory can be used at a time.