#### MODEL

Jazzy Elite ES-2S-SS	K0822
Jazzy Elite ES-2S-C	К0823
Weight Capacity	300 lbs.
Turning Radius	24.75″
Overall Width <sup>2</sup>	22.75″
Overall Length <sup>2</sup>	34.75″
Maximum Speed <sup>1</sup>	Up to 4 mph
Ground Clearance <sup>2</sup>	1.5″
Front Wheels	3" solid anti-tips
Drive Wheels	9" solid
Rear Wheels	6" solid casters
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Suspension Type	Limited
High-back Seat Sizes	16' x 16"-18", 18" x 18"-20", 20" x 18"-20"
Seat-to-Floor Range	21.5" - 23.5"
Per Charge Range <sup>2,4</sup>	Up to 13.5 miles
Standard Electronics	40A, PG GC3
Component Weights	Base 69 lbs. <i>(heaviest piece)</i> Seat 43 lbs. (comfort high-back) Batteries 25 lbs. each
Battery Requirements <sup>3,5</sup>	(2) 12 volt, U-1 batteries
Battery Charger	3.5 amp, off-board
Warranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

- 2 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).
- 3 AGM or gel cell type required.
- 4 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

5 Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or -3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



## **OPTIMAL** INDOOR MANEUVERABILITY IN SMALL SPACES



AMERICA'S FAVORITE POWER CHAIR gojazzy.com | 🖸 🔂 🕄



### (US) 800-800-8586 • (Canada) 888-570-1113 pridemobility.com

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice. JAZZYELITEES\_112119

# POWER, PERFORMANCE AND STYLE

Depth-adjustable, limited-recline, fold-flat, high-back seat

Intuitive controller is easy to operate

Flip-up foot platform offers greater comfort



Front-wheel drive design






PG GC3 Controller

An experience so superior, it can only be elite. The Jazzy<sup>®</sup> Elite ES features large, front-drive wheels to absorb the bumps and uneven surfaces of everyday life. The comfortable, high-back seat gives you the support you need all day, every day.

- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- Easy-access free-wheel lever
- 60" lap belt

#### **Available Colors**



Viper Blue

Candv Apple Red

#### **Optional Seating**

Deluxe contour, limited recline, high back with solid seat pan, available sizes:

- 16"W x 16"D
- 16"W x 18"D
- 18"W x 18"D
- 18"W x 20"D

Deluxe contour, limited recline, available sizes:

• 18"W x 16"D

Deluxe contour, limited recline, high back with solid seat pan, available sizes:

• 16"W x 16"D or 18"W x 20"D

Depth-adjustable, comfort seat, limited recline highback with solid seat pan, available sizes:

• 16"W x 16" - 18"D or 18"W x 18" - 20"D

#### Options

- Bell-shaped foot platform
- Elevating leg rests
- Heel loops
- Height-adjustable foot platform bracket
- Residual limb support
- Swing-away legrests
- Swing-away joystick

#### Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger