

The Traveler L4 is a high strength, lightweight Medicare approved K0004 wheelchair designed for both the home and rehab markets. Weighing less than 34 lbs, it features a dual position axle, multi-position caster forks for 18 or 20 inch seat height adjustment and offers flip-back desk and full length arm options for easy lateral movement. Precision bearings and non-marking polyurethane rear tires also expand the longevity of the wheelchair.



Height Adjustable Arms
(on select models)



Flip-Back Arms

Standard Features

- Available in 16", 18", 20" seat widths. Also 18" and 20" widths with 18" seat depth
- Quick release axle is available for some models
- Weighs less than 34 lb without front rigging
- Hammertone paint provides a long lasting and durable finish
- Carbon steel frame provides strength and durability
- Black nylon upholstery with standard chart pocket on back
- Anti-tippers come standard
- Tool-free adjustable front rigging
- Non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Flip-back arms
- Height adjustable flip-back arms on selected models
- Height adjustable back
- Two position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb (300 lb on models with 18" seat depth) - EVENLY DISTRIBUTED
- Meets Medicare code K0004



Warranty

Five (5) year limited warranty

**This product is not made with natural rubber latex*



Chart Pocket



Multi-position Caster













Adjustable Front Rigging



Industry Compatible Front Rigging



Anti-tippers come standard

Traveler® L4		Dimensions and Weights									
		Seat			Back	Overall				Wt.	Max Wt.
											
3F020120	Flip-Back Desk Arm, Swingaway Footrest	18"	16"	18" or 20"	16½" or 18½"	31½"	36"	25"	11½"	31.8 lb	250 lb
3F020130	Flip-Back Desk Arm, Elevating Legrest	18"	16"	18" or 20"	16½" or 18½"	31½"	36"	25"	11½"	31.8 lb	250 lb
3F020140	Flip-Back Full Arm, Swingaway Footrest	18"	16"	18" or 20"	16½" or 18½"	31½"	36"	25"	11½"	31.8 lb	250 lb
3F020150	Flip-Back Full Arm, Elevating Legrest	18"	16"	18" or 20"	16½" or 18½"	31½"	36"	25"	11½"	31.8 lb	250 lb
3F020220	Flip-Back Desk Arm, Swingaway Footrest	16"	16"	18" or 20"	16½" or 18½"	31½"	36"	23"	11½"	31.6 lb	250 lb
3F020230	Flip-Back Desk Arm, Elevating Legrest	16"	16"	18" or 20"	16½" or 18½"	31½"	36"	23"	11½"	31.6 lb	250 lb
3F020320	Flip-Back Desk Arm, Swingaway Footrest	20"	16"	18" or 20"	16½" or 18½"	31½"	36"	27"	11½"	32 lb	250 lb
3F020330	Flip-Back Desk Arm, Elevating Legrest	20"	16"	18" or 20"	16½" or 18½"	31½"	36"	27"	11½"	32 lb	250 lb
3F020340	Flip-Back Full Arm, Swingaway Footrest	20"	16"	18" or 20"	16½" or 18½"	31½"	36"	27"	11½"	32 lb	250 lb
3F020350	Flip-Back Full Arm, Elevating Legrest	20"	16"	18" or 20"	16½" or 18½"	31½"	36"	27"	11½"	32 lb	250 lb
3F030120	Flip-Back Desk Arm, Swingaway Footrest	18"	18"	18" or 20"	16½" or 18½"	33"	36"	25"	11½"	32 lb	300 lb
3F030130	Flip-Back Desk Arm, Elevating Legrest	18"	18"	18" or 20"	16½" or 18½"	33"	36"	25"	11½"	32 lb	300 lb
3F030140	Flip-Back Full Arm, Swingaway Footrest	18"	18"	18" or 20"	16½" or 18½"	33"	36"	25"	11½"	32 lb	300 lb
3F030320	Flip-Back Desk Arm, Swingaway Footrest	20"	18"	18" or 20"	16½" or 18½"	33"	36"	27"	11½"	32.2 lb	300 lb
3F030330	Flip-Back Desk Arm, Elevating Legrest	20"	18"	18" or 20"	16½" or 18½"	33"	36"	27"	11½"	32.2 lb	300 lb
With Quick Release Wheels											
3F020160	Adj. Height Desk Arm, Swingaway Footrest	18"	16"	18" or 20"	16½" or 18½"	31½"	36"	25"	11½"	31.8 lb	250 lb
3F020170	Adj. Height Desk Arm, Elevating Legrest	18"	16"	18" or 20"	16½" or 18½"	31½"	36"	25"	11½"	31.8 lb	250 lb
3F020260	Adj. Height Desk Arm, Swingaway Footrest	16"	16"	18" or 20"	16½" or 18½"	31½"	36"	23"	11½"	31.6 lb	250 lb
3F020270	Adj. Height Desk Arm, Elevating Legrest	16"	16"	18" or 20"	16½" or 18½"	31½"	36"	23"	11½"	31.6 lb	250 lb
3F020360	Adj. Height Desk Arm, Swingaway Footrest	20"	16"	18" or 20"	16½" or 18½"	31½"	36"	27"	11½"	32 lb	250 lb
3F020370	Adj. Height Desk Arm, Elevating Legrest	20"	16"	18" or 20"	16½" or 18½"	31½"	36"	27"	11½"	32 lb	250 lb
3F030160	Adj. Height Desk Arm, Swingaway Footrest	18"	18"	18" or 20"	16½" or 18½"	31½"	36"	25"	11½"	32 lb	300 lb
3F030170	Adj. Height Desk Arm, Elevating Legrest	18"	18"	18" or 20"	16½" or 18½"	31½"	36"	25"	11½"	32 lb	300 lb
3F030360	Flip-Back Desk Arm, Height Adj. Arm, Swingaway Footrest	20"	18"	18" or 20"	16½" or 18½"	33"	36"	27"	11½"	32.2 lb	300 lb

Traveler® L4 – Common Accessories & Parts		
90763064	Heel loop, pair	E0951
90763070	Push to lock, pair (includes rubber tips)	
90763071	Pull to lock, pair (includes rubber tips)	
90763062	Wheel lock extension, pair	E0961
300-505	Positioning belt (auto-style) push button	E0978
92133040	Positioning belt (airline style) (metal buckle)	E0978
B10320F	O2-IV holder	E2208
B10325	O2 cylinder holder	E2208
EJO2BAG07	O2 holder bag	E2208
90763230	Anti-tippers adjustable, pair	E0971
92523850	Anti-fold bar (chromed) lower frame	
9B761770	Anti-theft cap for legrest for Traveler chairs with 1.5" spacing	